


# Instrument Tables

**IT 2S, IT 3F, IT 3L, IT 4L,  
IT SPT, IT V.**





When your patients entrust you with their eyesight, their vision and your expertise converge. For the world's most advanced surgical and diagnostic solutions in ophthalmology, you can turn to Carl Zeiss Meditec. We're committed to earning your trust anew, every day.

# Technical Data.

Instrument tables	IT 2S	IT 3F	IT 3L	IT 4L	IT SPT	IT V
Footprint (max.)	700 x 480 mm	900 x 665 mm	950 x 650 mm	1000 x 650 mm	720 x 560 mm	900 x 665 mm
Table top	705 x 420 mm	900 x 420 mm	900 x 420 mm	1050 x 420 mm	620 x 550 mm	978 x 550 mm
Minimum table height	approx. 625 mm	approx. 625 mm	approx. 625 mm	approx. 625 mm	approx. 650 mm	approx. 667 mm
Vertical travel range	290 mm	290 mm	290 mm	290 mm	290 mm	290 mm
Vertical travel speed	17 mm/s	17 mm/s	17 mm/s	17 mm/s	17 mm/s	17 mm/s
Weight	19 kg	33 kg	45 kg	35 kg	28 kg	35 kg
Permissible load	15 kg	60 kg	20 kg	50 kg	40 kg	100 kg (centered on lifting column)
Rated voltage 120 V AC ± 10% 50 / 60 Hz	✓	✓	✓	✓	–	–
Rated voltage 230 V AC ± 10% 50 / 60 Hz	✓	✓	✓	✓	✓	✓

- Motorized height adjustment
- Outstanding stability
- Ergonomic design
- Suitable for patients in wheelchairs



**Carl Zeiss Meditec AG**  
 Goeschwitzer Str. 51-52  
 07745 Jena  
 Germany  
 Phone: +49 (0) 36 41 / 22 0-3 33  
 Fax: +49 (0) 36 41 / 22 0-2 82  
 info@meditec.zeiss.com  
 www.meditec.zeiss.com

**Carl Zeiss Meditec Inc.**  
 5160 Hacienda Drive  
 Dublin, CA 94568  
 USA  
 Toll Free: (800) 3 42-98 21  
 Phone: +1 (925) 5 57-46 51  
 Fax: +1 (925) 5 57-42 17  
 info@meditec.zeiss.com  
 www.meditec.zeiss.com