

DigiCam Adapter.
Digital, simple, versatile.





Your link to higher image quality.

The DigiCam adapter allows various digital cameras to be connected to Carl Zeiss Meditec slit lamps. Excellent images and unequalled operating comfort set a new standard in the field of slit lamp photography.

Ergonomic design:
Everything in view and under control

Positioning the camera directly beside the microscope simplifies operation. The operator hardly needs to move their head to have the LCD screen of the camera immediately in view. The camera controls are within easy reach, allowing camera setting changes quickly and conveniently, even during the examination.

Application:
Diverse options

Sony cameras with high-grade Carl Zeiss Vario Sonnar® objective and remote cable control are best suited for the DigiCam adapter. Therefore, the corresponding cameras, Sony DSC-S85 and Sony MVC-CD400, are included in our product range. For video documentation, the Sony DCR-PC120 digital camcorder is available.

The open architecture of the DigiCam adapter allows further digital cameras to be connected, such as the cameras of the popular Nikon CoolPix® Series.

Image quality:
Increased

Because of their high light transmission, Carl Zeiss Meditec slit lamps are optimally suitable for the documentation of findings. The DigiCam adapter, too, is equipped with high-grade Carl Zeiss optics. This results in brilliant pictures, whether survey, slit or high-power endothelium pictures.



**Additional illumination:
Bright without flash**

In taking slit images, the DigiCam illuminator ensures homogeneous illumination of the surrounding slit field, with the brightness being continuously adjustable. The result: First-class pictures without flash! Due to the use of a white-light LED, the reflection at the cornea is kept insignificant. The DigiCam illuminator does not need an additional power supply.

**Practical advantage:
Ready in no time**

The DigiCam adapter mounts to a standard beam splitter. This way, the SL 120 or the SL 130 diagnostic slit lamp easily turns into a digital photo slit lamp.

You can connect and disconnect the digital camera in no time and thus use it flexibly at different sites and for varied purposes.

Extension possibility: VISUPAC 121

As ideal complement to camera connection via the DigiCam adapter, we recommend the VISUPAC 121 digital image archiving and editing system. The VISUPAC software specially developed for ophthalmology includes a professional patient database and many image editing functions.





When your patients entrust you with their eyesight, their vision and your expertise converge. For the world's most advanced surgical and diagnostic solutions in ophthalmology, you can turn to Carl Zeiss Meditec. We're committed to earning your trust anew, every day.



Carl Zeiss Meditec AG
Goeschwitzer Str. 51-52
07745 Jena
Germany

Phone: +49 (0) 36 41 / 22 0-3 33
Fax: +49 (0) 36 41 / 22 0-2 82
info@meditec.zeiss.com
www.meditec.zeiss.com

Carl Zeiss Meditec Inc.
5160 Hacienda Drive
Dublin, CA 94568
USA

Toll Free: (800) 3 42-98 21
Phone: +1 (925) 5 57-46 51
Fax: +1 (925) 5 57-42 17
info@meditec.zeiss.com
www.meditec.zeiss.com