



FORUM Archive

The connection to a better communication - moving forward

The challenge

Many practices need an economic solution for sharing the diagnostic exams and patient demographic data between individual instruments as well as for safely storing and archiving these data. Diagnostic data needs to be available for review independently of the instrument.

The FORUM® solution

The FORUM Archive connects the diagnostic instruments and centrally stores their data, allowing you to share exams and patient demographic data. Diagnostic instruments from Carl Zeiss are also capable of storing their raw data in the FORUM Archive. This opens up the possibility e.g. for progression analyses across individual Field Analyzers (Humphrey HFA II-i) on one instrument. Furthermore, the FORUM Archive features a backup tool for data security.



The FORUM Archive connects your Carl Zeiss instruments and other DICOM compatible instruments and centrally stores their exams and patient demographic data.

Outstanding features

- Accessibility of the centrally stored data on all connected instruments
- Consistent patient data over all connected instruments
- One-time entry of patient demographic data to minimize mistakes and save time
- Remote analysis and review of diagnostic data via review software – independently from instruments
- Backup functionality providing data safety and integrity

Advanced version






Technical data

FORUM® Archive technical data. Minimum requirements:

CPU	<i>Intel Pentium 1.8 GHz (or more)</i>
RAM	<i>min. 2 GB</i>
Free hard disc capacity	<i>min. 3 GB</i>
Display resolution	<i>min. 1280 x 800 pixel</i>
Operating system	<i>Windows XP (Service Pack 3) or Windows Server 2003 (Service Pack 2) (32 bit) Windows Server 2008 R2 (64 bit) Windows 7 (64 bit)</i>

FORUM




 Carl Zeiss Surgical GmbH
 Carl-Zeiss-Strasse 22
 73447 Oberkochen, Germany

For more information, please visit:
www.meditec.zeiss.com/forum

Carl Zeiss Meditec AG
 Goeschwitzer Str. 51–52
 07745 Jena
 GERMANY

Phone: +49 36 41 22 03 33
 Fax: +49 36 41 22 01 12
info@meditec.zeiss.com
www.meditec.zeiss.com

Carl Zeiss Meditec Inc.
 5160 Hacienda Drive
 Dublin, CA 94568
 USA

Phone: +1 925 557 41 00
 Fax: +1 925 557 41 01
info@meditec.zeiss.com
www.meditec.zeiss.com