




FORUM

Configuration overview



FORUM is a modular system offering scalability and flexibility to fit the needs of your ophthalmic practice or clinic

The FORUM® Archive is the heart of all FORUM solutions. It connects ophthalmic instruments and centrally stores patient demographic data, diagnostic images, reports and raw data as well as surgery videos and reports. It supports the DICOM standard, thus offering compatibility with other DICOM instruments and systems. The FORUM Archive can also be linked to your practice management system (PMS) or electronic medical record (EMR) and allows

instrument-independent access of diagnostic data from Carl Zeiss review software applications.

The FORUM Viewer offers remote access and display of both static and dynamic data – allowing you to instantly view and flexibly compare all relevant data. Finally the new FORUM OR Manager provides surgery planning tools and relays all patient demographic and diagnostic data into the OR.

** CALLISTO eye connection to FORUM Archive requires FORUM OR Manager*

Advanced version





Technical data

FORUM® Archive technical data. Minimum requirements:

CPU	<i>Intel Pentium 1.8 GHz (or more)</i>
RAM	<i>min. 2 GB</i>
Free hard disc capacity	<i>min. 3 GB</i>
Display resolution	<i>min. 1280 x 800 pixel</i>
Operating system	<i>Windows XP (Service Pack 3) or Windows Server 2003 (Service Pack 2) (32 bit) Windows Server 2008 R2 (64 bit) Windows 7 (64 bit)</i>

FORUM Viewer technical data. Minimum requirements:

Hardware:	Windows	Macintosh
CPU	<i>Intel Pentium 1.8 GHz (or more)</i>	<i>1.8 GHz (or more)</i>
RAM	<i>min. 2 GB</i>	<i>min. 2 GB</i>
Free hard disc capacity	<i>min. 3 GB</i>	<i>min. 1 GB</i>
Display resolution	<i>min. 1280 x 800 pixel</i>	<i>min. 1280 x 800 pixel</i>
Operating system	<i>Windows XP (Service Pack 3) or Windows 7 (64 bit)</i>	<i>Snow Leopard (10.6.2)</i>

FORUM OR Manager technical data. Minimum requirements:

Hardware:	Windows	Macintosh
CPU	<i>Intel Pentium 1.8 GHz (or more)</i>	<i>1.8 GHz (or more)</i>
RAM	<i>min. 2 GB</i>	<i>min. 2 GB</i>
Free hard disc capacity	<i>min. 3 GB</i>	<i>min. 1 GB</i>
Display resolution	<i>min. 1024 x 768 pixel</i>	<i>min. 1024 x 768 pixel</i>
Operating system	<i>Windows XP (Service Pack 3)</i>	<i>Mac Leopard (10.5)</i>

FORUM



Carl Zeiss Surgical GmbH
Carl-Zeiss-Strasse 22
73447 Oberkochen, Germany

FORUM OR Manager



For more information, please visit:
www.meditec.zeiss.com/forum

Carl Zeiss Meditec AG
Goeschwitzer Str. 51–52
07745 Jena
GERMANY

Phone: +49 36 41 22 03 33
Fax: +49 36 41 22 01 12
info@meditec.zeiss.com
www.meditec.zeiss.com

Carl Zeiss Meditec Inc.
5160 Hacienda Drive
Dublin, CA 94568
USA

Phone: +1 925 557 41 00
Fax: +1 925 557 41 01
info@meditec.zeiss.com
www.meditec.zeiss.com