



FORUM EMR Gateway

Conveniently compatible with existing IT systems – the perfect fit

The challenge

Many eye care practices have an electronic medical record (EMR) system. To centrally store and access diagnostic data in the EMR system a connection between the EMR and the individual diagnostic instruments is required.

Our FORUM® solution

Using a simple EMR software upgrade, FORUM EMR Gateway creates a seamless link between

Carl Zeiss instruments and the existing EMR system.

This unique interface enables convenient transfer of diagnostic images and reports from the instrument to the EMR. Conversely, patient demographic data can be transferred from the EMR system to the instrument, including the CALLISTO eye OR Cockpit* from Carl Zeiss. To further facilitate efficient workflow, the optional FORUM Viewer can be launched from the EMR application, enabling automatic document selection and rapid, user-friendly display of current patient data (“desktop integration”).



Outstanding features

- Seamless link between Carl Zeiss instruments and the EMR system
- Easy transfer of instrument data to the EMR
- Transfer of patient demographic data from the EMR system to the instrument or to the CALLISTO eye OR Cockpit*
- One-time entry of patient demographic data to minimize mistakes
- Supports desktop integration when combined with FORUM Viewer
- Saves time and cost by eliminating the need for exam printouts

**CALLISTO eye connection to FORUM Archive requires FORUM OR Manager*





Technical data

FORUM® EMR Gateway technical data. Minimum requirements:


CPU	<i>Intel Pentium 1.8 GHz (or more)</i>
RAM	<i>min. 2 GB</i>
Free hard disc capacity	<i>min. 3 GB</i>
Display resolution	<i>min. 1280 x 800 pixel</i>
Operating system	<i>Windows XP (Service Pack 3) or Windows Server 2003 (Service Pack 2) (32 bit) Windows Server 2008 R2 (64 bit) Windows 7 (64 bit)</i>

FORUM EMR Gateway needs to be integrated by the EMR software provider.

Please contact your local ZEISS representative to learn which EMRs have already integrated FORUM EMR Gateway.

FORUM



 Carl Zeiss Surgical GmbH
Carl-Zeiss-Strasse 22
73447 Oberkochen, Germany

For more information, please visit:
www.meditec.zeiss.com/forum

Carl Zeiss Meditec AG
Goeschwitzer Str. 51-52
07745 Jena
GERMANY

Phone: +49 36 41 22 03 33
Fax: +49 36 41 22 01 12
info@meditec.zeiss.com
www.meditec.zeiss.com

Carl Zeiss Meditec Inc.
5160 Hacienda Drive
Dublin, CA 94568
USA

Phone: +1 925 557 41 00
Fax: +1 925 557 41 01
info@meditec.zeiss.com
www.meditec.zeiss.com