

# **ETHOS — First hydrophobic acrylic IOL injector**

ETHOS IOL Injector is the first injector preloaded with a hydrophobic implant injectable through a 2.8 mm incision. Held in a flat position, deflated and held by a safety device, the implant is packaged in the injector. This new injection system optimizes and facilitates the injection. Market launch in September 2006 at the ESCRS meeting in London.

- Safe: no handling guarantees perfect asepsis and the integrity of the implant
- Simple: the system is easy to use without extensive training
- Reproducible: the piston pushes the implant in a steady manner, standardising its exit from the injector into the eye
- Stable: the implant architecture is acknowledged for its performance in the capsular bag