



IHE Integration Statement			
Vendor	System Name	System Version	Date
Carl Zeiss Meditec	Cirrus HD-OCT	5.2	2011, Oct 10
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profile Implemented	Actor Implemented	Integration Profile Option Implemented	
Basic Eye Care Workflow	Acquisition Modality	None	
Internet address for vendors IHE information: http://www.meditec.zeiss.com/dcs			
Links to Standards Conformance Statements for the Implementation			
HL7			
DICOM	http://www.meditec.zeiss.com/dcs		
Links to general information on IHE			
In North America: www.ihe.net	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	