

SITA SWAP and HFA-Net Software Updates

Shipments of SITA SWAP and HFA-Net began in January 2006 for Humphrey® Field Analyzer customers who had these software products on order. Initial customer responses have been quite favorable.

Customers marvel that SITA SWAP runs so quickly, reducing testing time from 12-15 minutes per eye down to ~4 minutes per eye. Patients are much more willing to take a 4 minute test, and it fits easily into the clinical workflow.

HFA-Net brings a huge improvement for the electronic office by eliminating requirements for archiving and backing up to portable media, providing for seamless synchronization of HFA databases in the clinical practice, and printing to a network printer. A customer recently reported timesavings of ~2 minutes per eye per visual field after HFA-Net was installed. In a busy clinical practice this can have a significant impact on workflow efficiency.

NOTE: Installation of HFA-Net is a custom process. The task of establishing a successful network connection can be complex, and the need for a trained IT administrator is crucial. Customers must work with their onsite IT professional in order to install the software on their network.

Both SITA SWAP and HFA-Net are standard features on new HFA models and are also available as licensed upgrades for compatible HFA models. Please refer to the table below to view the software options available for the SITA SWAP and HFA-Net.

	SITA SWAP	HFA-Net
750i	x	x
745i	x	x
740i		x
720i		x