

VISULAS 532s and OPMI VISU 210 Surgical Microscopes Combine into One Unit

Both the OPMI VISU 210 surgical microscope and the VISULAS 532s laser offer a large number of innovative solutions for ophthalmic surgery. Now Carl Zeiss Meditec has combined these two devices into one space-saving unit. This was made possible by a new device support for the S88 floor stand. With it, the laser console can be set up on the operating microscope system.

The advantage for surgeons: a compact laser with exceptional performance requiring no extra space in the OR, and the cables to the operator protection filter are concealed in the stand arm. The control panel remains either on the console or is removed for flexible positioning. The operator also has more space with the 5m (16.4 feet) footswitch cable.

Carl Zeiss has leveraged its expertise in surgical and OR environments to develop the following:

- Improved ergonomics with Invertertube
- Fixed operator protection with minimal overall height

- Non-fatiguing, safe work environment with a laser display that is optimized for the work distance and the typical room lighting in the OR
- Sterile operation of the VISULAS 532s if necessary
- Surgery with a high-contrast, detailed and true-color microscope image with Superlux Eye xenon illumination

This powerful combination of OPMI VISU 210 and VISULAS 532s provide the optimal solution for better outcomes in the OR.

