



For the gynecologist's office: Colposcope E



Economy and quality in one.



Colposkop E - Low price, high quality -

The compact Colposcope E from Carl Zeiss proves that a flexible, very handy colposcope can also provide high optical quality.

And it also proves that "Economy" does not automatically mean less convenience. But it does mean an economical price!

The following models are available, all tailored precisely to your requirements in the gynecologist's office:

- on the suspension mount/on the floor stand
- with a straight tube/with an inclined tube
- with a fixed magnification/with a 3-step magnification changer
- video compatible:

The Medilive Compact Video Camera provides impressive information on the patient and superbly documents the results of the examination with the Colposkope E.

Technical Data

Illumination: 6V, 15W halogen illumination

Working distance: 300 mm


Magnification: 5x or with 3-step magnification changer: 3x - 5x - 8x

Green filter: swing-in

Suspension system: suspension mount (connecting piece required for examination chair*), mobile floor stand

*Mounts for attaching connecting pieces to examination chairs are available from the chair manufacturer.

Compliance:
DIN EN ISO 9001
EN 46001
ISO 13485



Subject to change in design and scope of delivery and as a result of ongoing technical development



*Application photos:
Prof. Hillemanns,
University Freiburg, Germany
Gynecological Hospital*

Carl Zeiss
Surgical Products Division
73446 Oberkochen
Germany
Fax: + 49 (0)7364/20-4823
email: surgical@zeiss.de
www.zeiss.de/colposcopy

