

EyeMag Smart - Medical Loupes
Seeing is believing.



EyeMag® Smart Medical Loupes
See the whole picture with stunning detail.

Easily recognizing structures helps you improve the quality of your treatment. EyeMag Smart loupes from Carl Zeiss provide optimum image quality at 2.5x magnification. The compact design of the optical components as well as the versatile adjustment possibilities promote comfortable working conditions and intuitive operation. EyeMag Smart loupes are a reliable partner for your day-to-day treatment.



EyeMag Smart

Optical systems					
Working distance (mm)	300	350	400	450	550
Working distance (inch)	12	14	16	18	22
Magnification	2.5x	2.5x	2.5x	2.5x	2.5x
Field of view (mm)	67	77	86	96	115

Compliance for EyeMag Smart Galilean Medical Loupes and accessories:



EyeMag is a registered trademark of Carl Zeiss.

Quality optics and sophisticated design

Precise vision...

- High image quality and color fidelity – brilliant image that extends to the peripheral zones
- High depth of field yields excellent depth orientation
- Large fields of view allow for an optimum overview of the working field

Intuitive operation ...

- Simple adjustment to individual requirements such as interpupillary distance and viewing angle
- Flip-up function facilitates unobstructed vision and easy communication with patients
- Easy placement and removal of loupe attachment streamlines daily use

Work in comfort ...

- Customizable settings enable you to sit in an ergonomically correct and comfortable position
- The large selection of working distances permits an optimal treatment position
- Even weight distribution minimizes fatigue during long periods of treatment

EyeMag Smart (Galilean type) Medical Loupes include:

- EyeMag Smart telescopes and carrier system (titanium eyeglass frame available in three sizes)
- Protective lens caps
- Sterilizable contact guards
- Optics cleaning cloth

Carl Zeiss Surgical GmbH
 A Carl Zeiss Meditec Company
 Carl-Zeiss-Straße 22
 73447 Oberkochen
 Germany
 Fax: +49 (0) 7364/20-4823
 Email: surgical@meditec.zeiss.com
www.meditec.zeiss.com/loupes
www.meditec.zeiss.com/contacts

