

**EyeMag Smart - Medical Loupes**  
**Seeing is believing.**



**EyeMag® Smart Medical Loupes**  
**See the whole picture with stunning detail.**

Easily recognizing structures helps you improve the quality of your treatment. EyeMag Smart loupes from Carl Zeiss provide optimum image quality at 2.5x magnification. The compact design of the optical components as well as the versatile adjustment possibilities promote comfortable working conditions and intuitive operation. EyeMag Smart loupes are a reliable partner for your day-to-day treatment.



*EyeMag Smart*

| Optical systems         |      |      |      |      |      |
|-------------------------|------|------|------|------|------|
| Working distance (mm)   | 300  | 350  | 400  | 450  | 550  |
| Working distance (inch) | 12   | 14   | 16   | 18   | 22   |
| Magnification           | 2.5x | 2.5x | 2.5x | 2.5x | 2.5x |
| Field of view (mm)      | 67   | 77   | 86   | 96   | 115  |

Compliance for EyeMag Smart Galilean Medical Loupes and accessories:

- DIN EN ISO 9001
- DIN EN ISO 13485
- CE

EyeMag is a registered trademark of Carl Zeiss.

**Quality optics and sophisticated design**

**Precise vision...**

- High image quality and color fidelity – brilliant image that extends to the peripheral zones
- High depth of field yields excellent depth orientation
- Large fields of view allow for an optimum overview of the working field

**Intuitive operation...**

- Simple adjustment to individual requirements such as interpupillary distance and viewing angle
- Flip-up function facilitates unobstructed vision and easy communication with patients
- Easy placement and removal of loupe attachment streamlines daily use

**Work in comfort...**

- Customizable settings enable you to sit in an ergonomically correct and comfortable position
- The large selection of working distances permits an optimal treatment position
- Even weight distribution minimizes fatigue during long periods of treatment

**EyeMag Smart (Galilean type) Medical Loupes include:**

- EyeMag Smart telescopes and carrier system (titanium eyeglass frame available in three sizes)
- Protective lens caps
- Sterilizable contact guards
- Optics cleaning cloth

**Carl Zeiss Surgical GmbH**  
 A Carl Zeiss Meditec Company  
 73446 Oberkochen  
 Germany  
 Fax: +49 (0) 7364/20-4823  
 Email: surgical@meditec.zeiss.com  
 www.meditec.zeiss.com/loupes  
 www.meditec.zeiss.com/contacts

