

EyeMag Light
Illuminate the Details.



EyeMag® Light

Powerful illumination for superior visualization.

The portable EyeMag Light is a powerful LED light source that brilliantly illuminates the visual field and ensures maximum detail recognition and high contrast. With its easy operation, integrated temperature regulation and a lightweight, compact design, EyeMag Light provides optimum wearing comfort for daily use, even during longer periods of treatment. The portable illumination system is an ideal supplement to the range of EyeMag loupes from Carl Zeiss.



Compliance for EyeMag® Light illumination source for loupes and accessories:

- DIN EN ISO 9001
- DIN EN ISO 13485
- DIN EN 60825-1
- LED class 2

EyeMag is a registered trademark of Carl Zeiss.

Carl Zeiss Surgical GmbH

A Carl Zeiss Meditec Company

73446 Oberkochen

Fax: +49 (0) 73 64 20-48 23

Email: surgical@meditec.zeiss.com

www.meditec.zeiss.com/loupes

www.meditec.zeiss.de/contacts



Powerful illumination and optimal operating comfort

Brilliant illumination for outstanding image quality

- Powerful LED light source ensures a bright illuminated field
- Optimum light distribution without color gradients
- Design and size of the illuminated field optimally adapted to EyeMag loupes
- Two brightness levels available

Simple operation and superior design for comfortable wear

- Easily attached to all EyeMag loupes
- Up to 8 hours of operating time
- Standard NiMH batteries for fast replacement
- Lightweight for maximum wearing comfort
- Portable battery supply unit with practical belt clip

Added safety features for reliable use during everyday treatment

- Medical device
- Active temperature regulation for optimal comfort
- LED indicator lights show the current operating status at all times
- Safe to sanitize unit surfaces with common disinfectants

EyeMag Light product line:

- EyeMag Light LED illumination
- Adapter for EyeMag Smart
- Adapter for EyeMag Pro

