

**CT ASPHINA™ 409M**  
**CT ASPHINA™ 409MV**



CT ASPHINA™ 409M



CT ASPHINA™ 409MV



## CT ASPHINA™ 409M CT ASPHINA™ 409MV



CT ASPHINA™ 409M



CT ASPHINA™ 409MV

Optic Design	Monofocal, aspheric (aberration neutral) Violet Light Filter (409MV)
Material	Hydrophilic acrylic (25%) with hydrophobic surface
Optic Diameter	6.0 mm
Total Diameter	11.0 mm
Haptic Angulation	0°
Lens Design	Single-Piece, MICS
Incision Size	1.5 mm
Company Labeled A-Constant*	118
Diopter Range	From 0.0 to +32.0 D  0.0 to +10.0 D, 1.0 D increments +10.0 to +30.0 D, 0.5 D increments +30.0 to +32.0 D, 1.0 D increments
ACD	4.97
Implantation in:	Bag
Injector / Cartridge Set**	AT.Shooter A2-2000 / ACM2 (1.5 mm)  Single-use injector A6 / *AT.Smart Cartridge Set (1.8 mm)

\* Please refer to our web pages for optimized A-Constants.

\*\* Please refer to our web pages for the most up-to-date references.

☞ Please ask your local representative for the current product name used in your country.

I/D01109 09/09 GB

ASPHINA™ is a trademark of Carl Zeiss Meditec.

The contents of the brochure may differ from the current status of approval of the product in your country.

Please contact our regional representative for more information.

Subject to change in design and scope of delivery and as a result of ongoing technical development.

© 2009 by Carl Zeiss Meditec SAS. All copyrights reserved.